



OPENCAST allows daily washing and showers with cast on. Water can be easily dried off with a towel and a hairdryer.



ODOR AND ITCH FREE

Perspiration is quickly

and unwanted odor.

evaporated to prevent itchiness

OPENCAST can be OPENCAST is made of fine repeatedly opened and be materials and has a stylish worn when needed. design.





S VISUAL CHECK

Doctors and patients can visually check the condition of the skin through windows while wearing the cast.



EASY REMOVAL

OPENCAST can be easily and safely removed.







open M inc.

413, Exporo, Yusung, Daejeon, 34051, Korea T 82.42.823.7301 | E sales@openm.com | www.openm.com






open M inc. OPENCAST OPENCAST is a new type of cast with a large net-like structure that allows

patients to take showers.

In addition, it is easy to fit and remove the cast, and patients will feel more comfortable and hygienic.

PRODUCT INFORMATION









CERVICAL COLLAR OPENCAST cervical collar can be customized to the patients' neck, resulting in several benefits for the customers.

BRACE

OPENCAST brace provides excellent air permeability. Due to the simple fitting process, brace can be easily made at the hospital.



OPENCAST size

	(unit : mm)
CM-03	130 x 342
CM-04	147 x 391
CM-05	167 x 435
CM-06	183 x 479
CM-07	204 x 527
CM-08	228 x 582

