

# Shower friendly cast



OPENCAST





# Instruction

**OPENCAST** is an ultimate water-friendly, customizable fixation solution. When heated, it can be formed into any anatomical shape, providing unlimited casting and customization possibilities. Its fully water-friendly design and big openings provide the patient and medical staff with an unrivaled level of convenience that will change the paradigm of casting once and for all.


## BENEFITS

- + Shower & Swim
- + No Itchiness & No odor
- + Breathable & Hygienic
- + Skin Monitoring




01

Heat **OPENCAST** using the oven to make it moldable. (soft sides facing each other)




03

Connect the two moldable plates of **OPENCAST** using clamps.




02

Put on stockinet of adequate thickness and width.




04

Mold **OPENCAST** around injured part and close the clamps.




05

Wrap with the transparent wrap, gently mold the shape, and let it cool down.



06

Unwrap the transparent wrap and remove the stockinet through the openings.



07

Check the final shape and for any patient inconveniences.



08

For foot, a **designated OPENCAST Shoe** must be always worn both indoor and outdoor except when washing or sleeping.

**openM inc.**

413, Expo-ro, Yuseong-gu, Daejeon, 34051, Korea

T +82.42.823.7301

E sales@openm.com

**www.openm.com**

v3.1 Nov. 2019